**parkell**<sup>®</sup>

EN

# Blu-Mousse<sup>®</sup> / Green-Mousse<sup>®</sup>

VPS Bite Registration / Impression Materials

## Instructions for Use

A00775rev10723

**REF** S457S, S448S, S440S, S438S, S437S, S436, S455S

Rx Only—product is intended for use by a licensed dental professional. For Safety Data Sheet (SDS) go to [www.parkell.com](http://www.parkell.com).

### Description:

Parkell's Blu-Mousse<sup>®</sup> and Green-Mousse<sup>®</sup> are addition-reaction silicones. (a.k.a. PVS, VPS.) When mixed, Mousse materials have a light, fluffy consistency, offer no resistance to closure, and have a rigid, hard final set. Set impressions are highly accurate, exhibit excellent tear strength, and have a dull surface appearance that is scannable without powder.

Mousse materials are available with set times of 30, 60 and 120 seconds. Dispensing options include dual barrel cartridges for gun dispensing, and tubes for hand mixing (*not available for all formulations*).

### Indications for Use:

Mousse materials are indicated for:

- Occlusal bite registrations
- Transfer impressions
- Instant bite rims
- Double-arch, single-stage or two-stage impressions (as the tray material).

### Storage and General Use Instructions:

- All Mousse products should be stored between 2 – 25°C (36 – 77°F). Exposure to excess heat will shorten shelf life and result in possible malfunction of the product.
- Avoid contact between Mousse materials and latex or sulfur-containing products, including latex gloves, rubber dams, and hemostatic agents with sulfates (e.g., Aluminum Sulfate, Ferric Subsulfate). If teeth or

gingiva are contaminated with these products, they must be thoroughly removed with copious amounts of water before impressing. Chloride containing hemostatic agents are acceptable substitutes.

- Vinyl polysiloxanes (VPS) are not compatible with polysulfides or polyethers. **DO NOT USE THEM TOGETHER.**
- Acrylic or methacrylate residues from bonding agents and restorative resins can inhibit impression setting. Completely remove this uncured layer from the surface of the prep and surrounding soft tissue prior to impressing.
- The set materials may be disinfected using glutaraldehyde or most other standard dental surface disinfectant products.
- Allergic reactions may occur in susceptible individuals. Flush affected tissues with copious water and consult a physician.
- Because of its high durometer, and plaster-like rigidity when set, to prevent impression lock-on, or the unintended removal of fixed restorations or teeth, NEVER use BluMousse to take impressions when the tooth and/or soft-tissue morphology exhibit the following characteristics: Undercuts, periodontally-involved teeth, pontics, and large embrasures.
- Never use rim-locking metal trays with Blu-Mousse. It will be extremely difficult to remove from the mouth if

the Blu-Mousse locks into undercuts, embrasures, or under pontics.

- Always evaluate and block out or lubricate all undercuts and fixed bridgework as appropriate before impressing, using wax, petroleum jelly or other VPS compatible materials.
- **IMPORTANT:** When using Blu-Mousse as the tray material for a double-step impression **DO NOT** use the Blu-Mousse impression as a matrix for bis-acryl temporary fabrication. The bis-acryl resin will contaminate the surface of the preliminary impression and prevent the second stage wash impression material from adhering to it
- **Use only genuine Parkell MIXPAC mixing tips, NOT counterfeit mixing tips!** Parkell is not responsible for damages caused by the use of unauthorized mixing tips.

## Using Gun Dispensed Mousse Materials:

1. Insert Mousse cartridge into the appropriate cartridge compatible impression gun. Remove the shipping cap but do not discard as you will be using it later to recap and reseal the cartridge.
2. Bleed about ¼" of material from both orifices, ensuring that it extrudes evenly from both openings
3. Attach a compatible mixer to the cartridge and twist-tighten until locked into position.

### **Tip/Cartridge Compatibility Guidelines:**

- Pink base for Split Cartridge (REF S445)
4. Clean rinse and dry the area being impressioned to remove all residual traces of any contaminants that may inhibit the setting of the material.
  5. Bleed a small amount of Mousse thru the mixer onto a pad then immediately dispense the material for the desired application.
  6. Wait the appropriate time (30, 60, or 120 seconds depending on the Mousse version used) and remove the set impression with a quick, firm motion.
  7. Do not reuse the shipping cap. Leave the used mixing tip on the cartridge as the new sealing cap.
  8. Mixers that were used intraorally should be removed from the cartridge and discarded. These cartridges should be resealed with the original sealing cap. To prevent cross contamination of the base and catalyst, the original cap should be cleaned of residual material before resealing. For mixers that were not used intraorally, the mixer can be left on the cartridge as a sealing cap.
  9. Parkell's Quick-Pour formula allows Mousse impressions to be poured in any dental stone.

### **Using Blu-Mousse Classic in Tubes:**

1. Dispense equal lengths of base and catalyst (1:1) onto a mixing pad. DO NOT Alter the mix ratio. It will impair final set.
2. Recap tubes immediately, making sure not to inter-change the color-coded caps.
3. Mix by hand with clean spatula until color is uniform and streak-free, with no air bubbles.
4. Load the material into an appropriate dispensing syringe for use.

### **Bite Registration Procedure with Mousse Materials:**

1. Mousse material may be dispensed directly from a mixing tip or syringe onto the biting surfaces of teeth, or into a conventional disposable plastic tray, side-less tray, or a "triple tray."
2. Always use tray adhesive when appropriate.

### **Tips for Dental Practitioners for Cross-Contamination Control:**

The following information should help dental practitioners reduce the risk of cross-contamination between patients when using multiple-use dental dispensers: apply disposable barrier sleeves/wraps over multiple-use dental dispensers before use with

each patient; use new, uncontaminated gloves when handling multiple-use dental dispensers; utilize dental assistants to dispense material for the dentist when appropriate; avoid contact of the reusable parts (e.g., the body of the multiple-use dental dispenser) with the patient's mouth; do not reuse the multiple-use dental dispenser if it becomes contaminated; do not immerse multiple use dental dispensers in a high level chemical disinfectant, as this may damage the dispenser and the material contained in the device; do not sterilize multiple-use dental dispensers, as this may damage the material contained in the device.

The cartridge and dispensing gun can be cleaned of debris with an EPA-registered low-level to intermediate-level hospital disinfectant. The used mixing tip should be discarded, and the original cap should be reattached to the cartridge to re-seal until next use. Utilize the surface disinfecting protocol of the disinfectant manufacturer. Consult [www.CDC.gov](http://www.CDC.gov) for the most recent version of the "Guidelines for Infection Control in Dental Health-Care Settings".

## Hardness Specification:

- **Blu-Mousse Hardness:**  
88 Durometer
- **Green Mouse Hardness:**  
60 Durometer

## Setting Time:

- **30 Seconds:** [REF](#) S457S, S448S
- **60 Seconds:** [REF](#) S440S
- **120 Seconds:** [REF](#) S438S, S436, S437S, S455S

## Disposal of Material:

Before disposal, the product should be fully mixed with the mixing tip, cured and then discarded in the general refuse. All cartridges and mixing tips may be disposed of in normal refuse.

## Events:

Reportable incidences should be reported to Parkell and the competent authority in which the user and or patient is established.

## Warranty and Terms of Use:

For full Warranty and Terms of Use information, please go to [www.parkell.com](http://www.parkell.com). Safety Data Sheets (SDS) are available at [www.parkell.com](http://www.parkell.com). Parkell's Quality System is certified to ISO 13485.

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## Explanation of Symbols Used:

|  |   |
|--|---|
|  | Professional use only                                   |
|  | Follow instructions for use                             |
|  | Temperature limit                                       |
|  | Store product tightly sealed away from direct sunlight. |
|  | Use-by date (expiration)                                |
|  | Batch code  |
|  | Catalogue / stock number                                |
|  | Unique Device Identifier                                |
|  | Medical Device  |

|   |                                  |
|---|----------------------------------|
|  | Package contents                 |
|  | Importer                         |
|  | Do not use if package is damaged |
|  | Single use only                  |
|  | Manufacturer                     |
|  | Translation                      |
|  | 50 ml Mixing Tip                 |
|  | 50ml Cartridge / Paste           |
|  | Base/Catalyst Tube / Paste       |

 **Parkell, Inc.**, 300 Executive Dr., Edgewood, NY 11717 USA  
(800) 243-7446 • [www.parkell.com](http://www.parkell.com)

**UK CA** **UKRP: Topdental (Products) Ltd,**  
12 Ryefield Way, Silsden, West Yorkshire BD20, 0EF England UK

 **mt-g medical translation GmbH & Co KG**  
Stuttgarter Straße 155 • D-89075 Ulm

  **MADE IN USA**