Panasil®

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Panasil[®]monophase medium

Panasil[®]initial contact regular

Panasil[®]initial contact regular fast

Panasil[®]contact two in one light

Panasil[®]initial contact light

Panasil[®]initial contact light fast

Panasil[®]contact plus x-light

Panasil[®]initial contact x-light

Panasil[®]initial contact x-light fast

Manufacturer:

Kettenbach GmbH & Co. KG

Im Heerfeld 7

35713 Eschenburg, Germany

www.kettenbach.com

Distributed in the U.S. by:

Kettenbach LP

7777 Center Avenue, Suite 280

Huntington Beach, CA 92647, USA

Vinyl Polysiloxane Impression Material ISO 4823

Made in Germany
37286/0810

For professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

Vinyl Polysiloxane Impression Material ISO 4823

Product Description

Panasil® impression materials are addition-curing, elastomeric materials with hydrophilic properties, high tear strength, dimensional accuracy, and a high resistance to permanent deformation. Panasil® initial contact displays a high initial hydrophilicity.

The Panasil® family consists of five different viscosities (putty, heavy-bodied, medium-bodied, light-bodied, x-light-bodied), available in an assortment of delivery systems: standard 1:1 50 ml automix cartridges, 5:1 362 ml foil bags for use in most automatic dispensing and mixing systems, and traditional 1:1 putty jars.

Indications/Techniques

Panasil® putty (putty, putty fast set, putty soft) and Panasil® binetics putty (putty fast, putty soft) are to be used as preliminary materials for:

- Two-step putty-wash impression technique
- One-step putty-wash impression technique

- Two-step putty-wash impression technique using a foil (plastic putty spacer)
- One-step putty impression technique for forming functional peripheries

Panasil® tray fast, Panasil® tray soft and Panasil® tray soft fast are to be used as heavy-bodied materials for:

- One-step impression technique (simultaneous technique) using single or dual viscosities
- Two-step impression technique using dual viscosities
- Functional impressions

Panasil® monophasic is to be used as a medium-bodied tray or syringeable impression material for:

- Taking impressions over fixed/removable restorations and implants (i.e., transferring impression posts and bridge components)
- Functional impressions
- Fabricating crown and bridgework or inlays

- Fabricating full or partial dentures
- Reline impressions
- Use in the simultaneous mixing technique as well as the putty-wash and triple tray techniques
- Transferring root posts when fabricating posts and cores indirectly

Panasil® contact plus, Panasil® contact two in one and Panasil® initial contact are to be used as syringeable impression materials for:

- Two-step putty-wash impression technique
- One-step putty-wash impression technique
- One-step impression technique using a foil (plastic putty spacer)
- One-step impression technique (simultaneous technique) using dual viscosities
- Reline impressions
- Fabricating full or partial dentures

Warnings

Do not use Panasil® impression materials as a temporary reliner.

Do not use Panasil® impression materials with condensation-curing silicones, polyether, Vinylsiloxan-ether® or polysulfide materials.

Panasil® putty impression materials are not suitable for detailed impressions when used alone.

Cautions

Do not use after expiration date.

Do not leave any residual impression material in the sulcus or oral cavity.

Do not swallow impression material! If swallowed: consult a medical doctor if problems arise or persist.

Avoid contact with eyes. If accidental contact with the eyes occurs, rinse immediately and thoroughly with an eye wash or water. Consult a medical doctor if problems arise or persist.

The product does not normally cause allergic reactions. However, for sensitive persons, an allergic reaction to the product cannot be ruled out.

Use of products containing active sulfur, aluminium chloride or nitrogen compounds (retraction cords containing ferric sulfate, polysulfide impression materials, etc.) in conjunction with this product will

interfere with the setting reaction of the vinyl polysiloxane material. Use of the materials requires the area to be rinsed thoroughly to remove all residue, before an impression is taken. Do not use latex gloves.

Do not interchange the base and hardener lids or scoops for hand-mixed jar putty.

In order to ensure an optimal adhesion of the two materials during the two-step putty-wash impression technique, both impression steps must take place one immediately after the other.

Before fitting the dynamic mixer, dispense material until equal amounts of base and catalyst appear; wipe excess. Firmly place dynamic mixer and lower locking lever.

When taking impressions of areas with severe undercuts and wide interdental spaces, use standard block-out techniques.

When using a custom impression tray, ensure that sufficient space remains between the side of the tray and the teeth/jaw.

Do not allow the material to enter the sewer or water system to avoid environmental contamination.

Avoid contact with clothing, since the material cannot be removed by dry cleaning.

Product Overview

Product name	ISO 4823	Consistency (approx.) mm	Mixing ratio and total content volume per unit	Mixing technique	Mixing element
Panasil® binetics putty fast	Type 0, Putty	23	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Panasil® binetics putty soft	Type 0, Putty	23	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Panasil® putty	Type 0, Putty	24	1:1, 900 ml in jars	manual, 1 scoop (12.3 ml) per component (Base & Hardener)	hand mix
Panasil® putty fast set	Type 0, Putty	24	1:1, 900 ml in jars	manual, 1 scoop (12.3 ml) per component (Base & Hardener)	hand mix
Panasil® putty soft	Type 0, Putty	24	1:1, 900 ml in jars	manual, 1 scoop (12.3 ml) per component (Base & Hardener)	hand mix
Panasil® tray fast heavy	Type 1, Heavy-bodied	32	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Panasil® tray fast heavy	Type 1, Heavy-bodied	34	1:1 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm
Panasil® tray soft heavy	Type 1, Heavy-bodied	32	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Panasil® tray soft heavy	Type 1, Heavy-bodied	34	1:1 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm
Panasil® tray soft heavy fast	Type 1, Heavy-bodied	32	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Panasil® tray soft heavy fast	Type 1, Heavy-bodied	34	1:1 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm

Product Overview

Product name	ISO 4823	Consistency (approx.) mm	Mixing ratio and total content volume per unit	Mixing technique	Mixing element
Panasil® monophase medium	Type 2, Medium-bodied	35	5:1, 362 ml in foil bags	Plug & Press® Dispenser, Sympress or other automatic dispensing and mixing unit	dynamic mixer
Panasil® monophase medium	Type 2, Medium-bodied	35	1:1 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm
Panasil® initial contact regular	Type 2, Medium-bodied	37	1:1 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm
Panasil® initial contact regular fast	Type 2, Medium-bodied	37	1:1 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm
Panasil® contact two in one light	Type 3, Light-bodied	38	1:1 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static green mixing tip MB Ø 6.5 mm
Panasil® initial contact light	Type 3, Light-bodied	41	1:1 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static yellow mixing tip MB Ø 4.2 mm
Panasil® initial contact light fast	Type 3, Light-bodied	41	1:1 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static yellow mixing tip MB Ø 4.2 mm
Panasil® contact plus x-light	Type 3, Light-bodied	42	1:1 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static yellow mixing tip MB Ø 4.2 mm
Panasil® initial contact x-light	Type 3, Light-bodied	44	1:1 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static yellow mixing tip MB Ø 4.2 mm
Panasil® initial contact x-light fast	Type 3, Light-bodied	44	1:1 50 ml cartridge	Applyfix® 4 dispensing gun DS-50 1:1/2:1	static yellow mixing tip MB Ø 4.2 mm

Technical Data

Product name	Mixing ratio	Working time at 23° C/ 74° F ≤	Working time at 35° C/ 95° F ≤	Intraoral setting time at 35° C/ 95° F ≥	Total setting time* ≥	Hardness (approx.) Shore	Linear dimensional change (approx.) %	Elastic-Recovery test (approx.) %	Strain in compression (approx.) %
Panasil® binetics putty fast	5:1	1 minute 30 seconds	not applicable	2 minutes 30 seconds	4 minutes	A 63	- 0.20	99.5	2.5
Panasil® binetics putty soft	5:1	2 minutes	not applicable	3 minutes	5 minutes	A 56	- 0.20	99.5	3.5
Panasil® putty	1:1	2 minutes	not applicable	2 minutes	4 minutes	A 66	- 0.20	99.0	2.7
Panasil® putty fast set	1:1	1 minute 30 seconds	not applicable	2 minutes	3 minutes 30 seconds	A 66	- 0.20	99.0	2.7
Panasil® putty soft	1:1	2 minutes	not applicable	2 minutes	4 minutes	A 60	- 0.20	99.0	2.7
Panasil® tray fast heavy	5:1	1 minute 20 seconds	not applicable	2 minutes	3 minutes 20 seconds	A 62	- 0.20	99.7	2.5
Panasil® tray fast heavy	1:1	1 minute	not applicable	2 minutes	3 minutes	A 62	- 0.20	99.5	3.0
Panasil® tray soft heavy	5:1	2 minutes	not applicable	2 minutes	4 minutes	A 55	- 0.20	99.7	3.0
Panasil® tray soft heavy	1:1	2 minutes	not applicable	2 minutes	4 minutes	A 55	- 0.20	99.5	3.0
Panasil® tray soft heavy fast	5:1	1 minute 20 seconds	not applicable	2 minutes	3 minutes 20 seconds	A 55	- 0.20	99.5	3.0
Panasil® tray soft heavy fast	1:1	1 minute	not applicable	2 minutes	3 minutes	A 55	- 0.20	99.5	3.0

*Total setting time (removal time from mouth) from start of mix.

Technical Data

Product name	Mixing ratio	Working time at 23° C/ 74° F ≤	Working time at 35° C/ 95° F ≤	Intraoral setting time at 35° C/ 95° F ≥	Total setting time* ≥	Hardness (approx.) Shore	Linear dimensional change (approx.) %	Elastic-Recovery test (approx.) %	Strain in compression (approx.) %
Panasil® monophase medium	5:1	2 minutes	1 minute	2 minutes	4 minutes	A 58	- 0.20	99.7	3.0
Panasil® monophase medium	1:1	2 minutes	1 minute	2 minutes	4 minutes	A 58	- 0.20	99.7	3.5
Panasil® initial contact regular	1:1	1 minute 30 seconds	1 minute	2 minutes 30 seconds	4 minutes	A 46	- 0.20	99.7	3.0
Panasil® initial contact regular fast	1:1	1 minute	30 seconds	2 minutes	3 minutes	A 46	- 0.20	99.7	3.0
Panasil® contact two in one light	1:1	2 minutes	1 minute	2 minutes	4 minutes	A 46	- 0.20	99.7	3.5
Panasil® initial contact light	1:1	1 minute 30 seconds	1 minute	2 minutes 30 seconds	4 minutes	A 46	- 0.20	99.3	3.5
Panasil® initial contact light fast	1:1	1 minute	30 seconds	2 minutes	3 minutes	A 46	- 0.20	99.3	3.5
Panasil® contact plus x-light	1:1	2 minutes	1 minute	2 minutes	4 minutes	A 46	- 0.20	99.7	5.0
Panasil® initial contact x-light	1:1	1 minute 30 seconds	1 minute	2 minutes 30 seconds	4 minutes	A 46	- 0.20	99.3	3.5
Panasil® initial contact x-light fast	1:1	1 minute	30 seconds	2 minutes	3 minutes	A 46	- 0.20	99.3	3.5

*Total setting time (removal time from mouth) from start of mix.

Note

To ensure optimal impressions, the temperature of the material should not deviate from 23° C (74° F) before applying. Otherwise, working and setting times will be affected.

Impression tray: Preparation and adhesives

In principle, all common impression trays can be used if a relevant dynamic pressure is guaranteed. When retention is not sufficient, ensure a strong bond to the impression material, by brushing the impression tray with a thin film of Panasil® adhesive for addition-curing silicone prior to loading the tray with impression material. Allow to dry per manufacturer's instructions.

Instructions for use: Jar material for hand-mix putty preparation

Only use materials with the same lot number.

Do not interchange the base and hardener lids or scoops.

Use the scoops to dispense equal amounts of base and hardener paste.

The different color scoops must only be used for the materials with corresponding colors.

Close the containers carefully after use, and ensure that the lids are not interchanged.

Contamination of base paste with hardener paste in the container renders the material unusable.

Knead the base paste and hardener paste for 45 seconds until the material is homogeneously colored.

Should gloves be worn, test them for compatibility prior to mixing with a sample of the material to be mixed. Certain types of gloves (such as latex) can prevent the polymerization. Using vinyl gloves is recommended.

Prior to taking the wash impression, the preliminary impression must be cleaned with water, dried and trimmed using standard methods.

Preparation of the Plug & Press® foil bag system

Foil bag material for use in most automatic dispensing and mixing systems

Before using the base material for the first time, remove the white safety pin from the activating head of the large foil bag by turning it in the direction of the arrow and pulling it out (Figure 1, Page 114/115).

Place the large foil bag with the integrated activating head in the cartridge body. Ensure that the notches on the activating head and cartridge case are aligned (Figure 2, Page 114/115).

Firmly press the activating head into the final position on the cartridge body. The foil bag is automatically pierced by a pin when the activating head is pressed down (Figure 3, Page 114/115).

Follow the same procedure as in Figure 2 and 3 for the catalyst material. (Note: There is no safety pin on the catalyst foil bag.)

The use of dynamic mixers and cartridge bodies supplied by Kettenbach is a prerequisite for an optimum result.

Application using the Plug & Press® Dispenser in combination with the Kettenbach's dynamic mixing tip

- Move the plungers to the top by turning the control knob (Figure 4, Page 114/115).
- Insert the cartridge body into the unit (Figure 5, Page 114/115).
- Turn the control knob to move the plungers into the cartridge body, and continue turning until the plungers come into contact with the foil bags (Figure 6, Page 114/115).
- Once there is contact, press one of the two start keys (at the front on top of the unit) to dispense the material (Figure 7, Page 114/115).
- Before fitting the dynamic mixer, dispense material until equal amounts of base and catalyst appear (Figure 7, Page 114/115); wipe excess. Firmly place dynamic mixer (Figure 8, Page 114/115) and lower locking lever.

Load the tray with the required amount of material. Hold the impression tray at an angle and press lightly against the tray. Leave the dynamic mixer in

the material while dispensing (Figure 9, Page 114/115). Leave the filled dynamic mixer on the cartridge body as a seal.

Before next use, release the locking clip to remove the used dynamic mixer and check that the outlets in the activating heads are not blocked. Fit a new dynamic mixer, lower locking lever and continue as usual.

When finished, simply remove the empty foil bags and activating heads from the cartridge body and discard them. The cartridge body is reusable (Figure 10, Page 114/115). The activating heads are disposable.

If not using a Kettenbach Plug & Press® Dispenser, please follow the manufacturer's instructions for your automatic dispensing and mixing system.

Application using the Sympress Dispenser in combination with the Kettenbach's dynamic mixers

- Allow the plungers (with cover closed) to return to the original position by pressing the plunger return button (Figure 1, Page 116/117). Next, open the cover (Figure 2, Page 116/117).

- Insert the cartridge body with the foil bags into the device (if foil bag material is already in use, the filled mixer will be on the cartridge body as the closure) (Figures 2 and 3, Page 116/117)
- Mount the dynamic mixing tip in place (if you are already using foil bag material, first remove the mixer, which is being used as the closure and replace it with the new dynamic mixing tip). (Figure 4, Page 116/117)
- Close the locking lever (Figure 5, Page 116/117)
- Close the device cover (Figure 6, Page 116/117)
- Dispense the material, dispose of the initial approx. 3 cm (until a homogenous mix of base and catalyst appear; see coloring of the final impressions material) (Figures 7 and 8, Page 116/117)
- Load the tray (or application syringe) (Figure 9, Page 116/117). As you do, hold the impression tray at an angle and press lightly against it. Leave the mixing tip in the material while dispensing to prevent air from getting trapped.

- Leave the filled dynamic mixer tip in place as a seal.

Before next use, release the locking clip to remove the used dynamic mixing tip and check that the outlets in the activating heads are not blocked. Fit a new dynamic mixing tip, lower locking lever and continue as usual.

When finished simply remove the empty foil bags and activating heads from the cartridge body and discard them. The cartridge body is reusable (Figure 10, Page 114/115). The activating heads are disposable.

If you are using an automated dosage and mixing system other than the Sympress Dispenser, please follow the instructions for your respective system.

Instructions for use:

Cartridge material for use in Applyfix® 4 dispensing gun DS-50 1:1/2:1

Insert the cartridge into the Applyfix® 4 dispensing gun DS-50 1:1/2:1. Ensure that the notches on the cartridge base are pointing down. The clasp will not close if the cartridge is not properly inserted into

the dispensing gun.

Remove the cartridge cap. The cap can be replaced after initial use.

Advance plunger into cartridge and dispense a small amount of impression material until equal amounts are extruded at the same rate.

Install a mixing tip on the cartridge, and turn cap 90 degrees clockwise to lock in place.

Inject the required amount of material, directly into the tray or the preliminary impression, or into a delivery syringe. To apply the material around the prepared teeth, use the delivery syringe or the cartridge with an intraoral tip on the mixer. Check to ensure materials are set prior to mouth removal.

Leave used mixing tip on the cartridge after use or replace the cartridge cap. Before using the cartridge again, remove cartridge cap or remove and discard the used mixing tip.

Check the cartridge openings to ensure no polymerized material is present. Should this be the case, unblock the openings and dispense a small amount of impression material until equal amounts are dispensed at the same rate.

Install a new mixing tip and continue as described above.

The Applyfix® 4 dispensing gun DS-50 1:1/2:1 can be sterilized in the autoclave.

Storage conditions:

Store in a dry place at room temperature; do not expose to sunlight and avoid exposure to extreme temperature fluctuations.

Disinfection:

Impressions can be disinfected using for example, a 2% acidic glutaraldehyde solution. Use disinfection solutions specific for impression materials (e.g., Silosept®). See manufacturer instructions.

Electroplating:

Impressions may be silver-plated with an alkaline-plating solution only.

Model preparation:

For model fabrication, the impression can be poured immediately after disinfection. Impressions can also be poured within two weeks with standard dental stone class IV (e.g., Tewestone® or Tewelock®).

Trademarks

Kettenbach®, Applyfix® 4, Panasil®, Silosept®, Tewelock®, Tewestone® and Plug & Press® Dispenser are registered trademarks of Kettenbach GmbH & Co. KG.

Note

Selected Kettenbach materials are available in certain markets only.

Limitation of liability

Except where prohibited by law, Kettenbach GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.



The information provided for Kettenbach products is based on comprehensive research and experience in application technology.

Results are furnished to the best of our knowledge, subject to technical changes within the framework of product development.

However, users must comply with and consider all recommendations and information in connection with any use.

Updated: 1st March 2011

Legend

Conformity with MDD 93/42/EEC		Use by	
Temperature limitation		Keep dry	
Keep away from sunlight		Catalog number	REF
Batch code		Caution, consult accompanying documents	
Dual system (Germany only)		Diameter	
Cartridge		Static mixing tip MB (single use)	
Mixing tip Type B	MB	Intraoral tip (single use)	
Scoop		Jar	
Adhesive		Single use	
Cartridge body		Foil bags	
Dynamic mixer (single use)		Rx only For professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.	
Millimeter	mm	Milliliter	ml

Panasil®

Panasil[®]binetics putty fast

Panasil[®]binetics putty soft

Panasil[®]putty

Panasil[®]putty fast set

Panasil[®]putty soft

Panasil[®]tray fast heavy

Panasil[®]tray soft heavy

Panasil[®]tray soft heavy fast

Panasil[®]monophase medium

Panasil[®]initial contact regular

Panasil[®]initial contact regular fast

Panasil[®]contact two in one light

Panasil[®]initial contact light

Panasil[®]initial contact light fast

Panasil[®]contact plus x-light

Panasil[®]initial contact x-light

Panasil[®]initial contact x-light fast

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Panasil®

Panasil[®]binetics putty fast
Panasil[®]binetics putty soft
Panasil[®]putty
Panasil[®]putty fast set
Panasil[®]putty soft
Panasil[®]tray fast heavy
Panasil[®]tray soft heavy
Panasil[®]tray soft heavy fast
Panasil[®]monophase medium
Panasil[®]initial contact regular
Panasil[®]initial contact regular fast
Panasil[®]contact two in one light
Panasil[®]initial contact light
Panasil[®]initial contact light fast
Panasil[®]contact plus x-light
Panasil[®]initial contact x-light
Panasil[®]initial contact x-light fast

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Panasil®

Panasil[®] binetics putty fast
Panasil[®] binetics putty soft
Panasil[®] putty
Panasil[®] putty fast set
Panasil[®] putty soft
Panasil[®] tray fast heavy
Panasil[®] tray soft heavy
Panasil[®] tray soft heavy fast
Panasil[®] monophase medium
Panasil[®] initial contact regular
Panasil[®] initial contact regular fast
Panasil[®] contact two in one light
Panasil[®] initial contact light
Panasil[®] initial contact light fast
Panasil[®] contact plus x-light
Panasil[®] initial contact x-light
Panasil[®] initial contact x-light fast

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