CONTENT

EN (English) Instruction Manual GenTek™ Pre-milled Abutment Blank.....	2
BG (Български) Ръководство за употреба Предварително фрезована заготовка за абатмънт GenTek™.....	5
CS (Česky) Návod k použití GenTek™ Předfrézovaná forma abutmentu.....	8
DA (DANISH) Brugervejledning GenTek™ forfræset abutmentemne.....	11
DE (Deutsch) Gebrauchsanweisung GenTek™ Pre-milled Abutment Blank.....	14
EL (Ελληνικά) Εγχειρίδιο οδηγιών Προκατασκευασμένο κολόβωμα GenTek™.....	17
ES (Español) Manual de instrucciones Blank para pilar prefresado GenTek™.....	20
ET (Eesti) Kasutusjuhend GenTek™-i i eelfreesitud abutmenditoorikud.....	23
FI (Suomi) Käyttöohje Esijyrsitty GenTek™ -abutmenttiaihio.....	26
FR (Français) Mode d'emploi GenTek™ Pre-milled Abutment Blank.....	29
HR (Hrvatski) Priručnik GenTek™ prethodno glodani valjčić za izradu nadogradnje.....	32
HU (Magyar) Kezelési útmutató GenTek™ Előmart Nyers Támpillér.....	35
IT (Italiano) Manuale d'istruzioni GenTek™ Pre-milled Abutment Blank.....	38
LT (Lietuvių k) Naudojimo instrukcija GenTek™ iš anksto išfrezuotas atramos ruošinys.....	41
LV (Latviski) Lietošanas rokasgrāmata GenTek™ iepriekš frēzēts abatments.....	44
NL (Nederlands) Gebruikershandleiding Voorgefreesde GenTek™ -opbouwbasisvorm.....	47
NO (Norsk) Instruksjonshåndbok GenTek™ forhåndsfrest abutmentform.....	50
PL (Polski) Instrukcja stosowania Pre-milled GenTek™.....	53
PT (Português) Manual de Instruções Blank para pilar prefresado GenTek™.....	56
RO (Română) Manual de instrucțiuni Bloc de bont protetic prefabricat GenTek™.....	59
SK (Slovak/Slovakian) Návod na použitie Predfrézované nevytvárané abutmenty GenTek™.....	62
SL (Slovensčina) Navodila za uporabo Predrezkani polizdelek za abutment GenTek™.....	65
SV (Svenska) Bruksanvisning GenTek™ förfräst distansämne.....	68
TR (Türkçe) Kullanım Kılavuzu GenTek™ Önceden Frezelenmiş Abutment Taslağı.....	71

Instruction Manual

GenTek™ Pre-milled Abutment Blank

including final prosthetic screw

Version: **3/06.2023**
 Compiled at: **23.06.2023**
 Effective date: **26.06.2023**
 Replaces version: **2/01.2022**

Indications

Intended Use / Intended Purpose

The GenTek™ Pre-milled Abutment Blank is intended to enable the production of a single-unit patient-specific abutment to support screw and cement-retained restoration.

Device description

Each GenTek™ Pre-milled Abutment Blank is packaged with the final prosthetic screw. The final prosthetic screw for each system is also available as an individually packaged component. The GenTek™ Pre-milled Abutment Blank and prosthetic screw are for single use only. The GenTek™ Pre-milled Abutment Blank and prosthetic screw are provided non-sterile and intended to be cleaned and sterilized prior to placement into the patient's mouth.

The GenTek™ Pre-milled Abutment Blank is available with genuine connections for the following implant-level Zimmer Biomet Dental implant systems: Certain®, External Hex, TSV®/Trabecular Metal™, and Eztetic®. For specific dental implant compatibility, please refer to the information on the individual product label. GenTek™ Pre-milled Abutment Blanks are compatible with Zimmer Biomet Dental prosthetic tooling as outlined in the GenTek™ Components Catalog – please refer to the components catalog located on the Zfx website.

Materials

The GenTek™ Pre-milled Abutment Blank is manufactured from biocompatible titanium alloy grade 5 ELI (TiAl6V4 ELI).

Procedure


× Designing the patient-specific abutment

The patient-specific abutment can be processed via a digital or traditional restorative workflow – please refer to the components catalog located on the Zfx website for compatible scan bodies. The patient-specific abutment must be designed using dental design software. The geometry of the pre-milled blank is contained in the library file. **Ensure that the connection of the pre-milled blank matches the connection and seating surface size of the implant.** Guidelines for recommended design parameters, minimum wall thickness and maximum angles are shown in Figure 1.

Guidelines for recommended design parameters:

For the design of the patient-specific abutment the following points should be considered:

1. The minimum wall thickness should be not less than 0.5 mm.
2. The maximum angulation (A) should not exceed 30 degrees.

× Milling a Custom Restoration

The machine must be able to process pre-milled abutment blanks. Confirm that the GenTek™ Pre-milled Abutment Blank library is installed in your CAM software (i.e. Hyperdent, Sum 3D, etc.). Select the appropriate blank as per the design. The blank must be properly assembled to the blank holder in the milling machine to ensure that the interface of the abutment blank is not damaged during the milling process. Please refer to www.zfx-dental.com for compatible CAM software and abutment blank holders.

Note: The position of the GenTek™ Pre-milled Abutment Blank must be calibrated using the GenTek™ calibration blank prior to machining. **The outcome of the milling process is a patient-specific abutment.**

× Placing / Removing the Patient-Specific Restoration

1. Remove any provisional abutment, cover screw, healing abutment or similar from the implant upon which the patient-specific abutment is intended to be placed.
2. Visually confirm that the selected patient-specific restoration is the same size as the implant.
3. Ensure that the seating surface of the implant is clean, free of residue, bone, or soft tissue and is dry when mounting.
4. Affix the patient-specific restoration into the implant with the prosthetic screw to the specified torque using the corresponding Zimmer Biomet screwdriver. Refer to the below table for final placement torque recommendations:

Connection	Final Prosthetic Screw Recommended Torque
Certain® (Gold-Tite® Screw)	20 Ncm
External Hex (Gold-Tite® Screw)	32 Ncm
TSV® / Trabecular Metal™	30 Ncm
Eztetic®	30 Ncm

Note: Unlike the Certain® and External Hex systems, the TSV® and Eztetic® Pre-Milled Abutment engaging connections form a friction-fit with the mating implant upon final assembly, requiring the use of additional instrumentation to remove the abutment from the implant.

× Removal of Friction Fit Connections

Should the patient-specific abutment (TSV®/Trabecular Metal™ and Eztetic®) restoration require disassembly from the implant upon loosening and removal of the screw utilizing a GemLock Tool, an abutment removal tool may be required. Refer to GenTek™ Friction-Fit Abutment Removal Tool IFU-4097 for abutment removal instructions.

Removal tool for narrow platform (NP – Eztetic)

Removal tool for standard platform (SP – TSV/TM)

Instruction Manual

GenTek™ Pre-milled Abutment Blank

including final prosthetic screw

Version: 3 / 06. 2023
Compiled at: 23. 06. 2023
Effective date: 26. 06. 2023
Replaces version: 2 / 01. 2022

Cleaning and Sterilization

The GenTek™ Pre-milled Abutment Blank is non-sterile. Prior to placement of the patient-specific abutment and prosthetic screw into the implant, the patient-specific abutment and prosthetic screw assembly must first be cleaned and sterilized.

Step 1: Manual Cleaning

1. Rinse the instruments under running cold tap water for 2 minutes.
2. Hand wash instruments in a bath with demineralized water for 5 minutes. Remove debris with a soft brush. Cavities should be rinsed with a syringe. This step should be repeated until all debris is removed.
3. Clean the instrument in an ultrasonic bath using a mildly alkaline cleaner (1.5%) at a frequency of at least 40 kHz for 20 minutes.
4. Rinse thoroughly with demineralized water for 1 min.

Step 2: Sterilization

One of the following sterilization methods is recommended to be used:

- × **Method 1:**
Full cycle pre-vacuum steam sterilization at a temperature of 134°C (273.2°F) for 3 minutes, 4 pulse minimum.
- × **Method 2: (alternative to method 1)**
Full cycle gravity steam sterilization at 134°C (273.2°F) for 3 minutes minimum.

Step 3: Drying

After sterilization treatment it is recommended to dry the products at a temperature of maximum 134°C (273.2°F) for 10 minutes. The products are ready for use 5 minutes after completion of the drying process.

Note: Follow the loading instructions issued by the sterilizer manufacturer.

Storage and Handling

The GenTek™ Pre-milled Abutment Blank should be stored at room temperature.

Known Side Effects

Allergies to the alloy or contents of the alloy may very rarely occur.

Warnings

- × The GenTek™ Pre-milled Abutment Blank must only be used and handled by dental professionals.
- × GenTek™ Pre-milled Abutment Blanks are not intended for use with non-Zimmer Biomet Dental implant systems.
- × The GenTek™ Pre-milled Abutment Blank and prosthetic screw must be handled with caution to avoid swallowing or inhaling by the patient.
- × Do not direct layer porcelain to the patient-specific abutment as this can affect mechanical properties of the titanium.
- × The use of the prosthetic screw torque higher than the recommended value may result in damage to the final abutment, the prosthetic screw and/or the implant. The use of torque lower than the recommended value may result in loosening of the final abutment and/or the prosthetic screw.
- × The screwdriver must not be worn or damaged. Use of the wrong or a worn or damaged screwdriver can lead to over-rotation of the hex in the screw head which can result in the screw not being fully tightened and/or not being able to be removed.
- × The included prosthetic screw may only be used to permanently fix the finished restoration. For lab use and fitting test a try-in screw must be used.

Instruction Manual

GenTek™ Pre-milled Abutment Blank

including final prosthetic screw

Version: 3 / 06. 2023
Compiled at: 23. 06. 2023
Effective date: 26. 06. 2023
Replaces version: 2 / 01. 2022

Explanation of symbols

Product is packaged non sterile

Do not use if package is damaged

Manufacturer Information: Address of legal manufacturer's facility

Don't use twice

Recommended torque to be applied for final placement of the included screw

Recommended to fix the included screw hand-tight

LOT-Number, Symbol followed by Production Date of the lot

Medical Device

Reference Number, Symbol followed by Manufacturer's Item Number

Quantity / Unit of package content

Unique Device Identifier

Refer to Instruction for use:
Document is available for download on Zfx-Homepage

Information on Manufacturer

Zfx GmbH +49 (0) 8131 / 33 244 - 0
Kopernikusstraße 15 F +49 (0) 8131 / 33 244 - 10
85221 Dachau www.zfx-dental.com
Germany

Incident Reporting / Complaints

Please send to emeacomplaints@ZimVie.com

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