Ecosite Elements

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English

Product description

Ecosite Elements is a light-curing, esthetic, radiopaque nanohybrid composite system for the anterior and posterior region. The range of highly viscous composites includes 1 bleach shade (BL1), 5 VITA shades (B1, A1, A2, A3 and A3.5) and 4 enamel shades (Enamel Bleach = EB, Enamel Light = EL, Enamel Medium = EM and Enamel Dark = ED) in different opacities. This range of shades, combined with the unique handling and excellent polishing properties of the composites, makes the creation of highly-esthetic restorations possible.

The shades A2O (A2 Opaque) and INC (Incisal) are radiopaque, flowable composites. The shade A2O is recommended for the veneering of slightly discolored areas and for preparing highly opaque teeth. The shade INC is particularly suitable for fashioning esthetic incisal edges thanks to its high transparency. Both shades are suitable for fissure sealing.

The shades W (White) and B (Brown) are flowable composites for the customization of composite restorations. The shade W is used for concealing white spots and as an opaque cover layer. The shade B is an ideal choice for fissures and the treatment of brown spots.

Intended purpose

Ecosite Elements, Shades: BL1, B1, A1, A2, A3, A3.5, Enamel Bleach, Enamel Light, Enamel Medium, Enamel Dark

- Restorations of all cavity classes of every dentition
- Direct veneers (vestibular and occlusal)
- Core build-up

Ecosite Elements, Shades: A2 Opaque and Incisal

- Restorations of the cavity Classes III, IV and V of every dentition
- Bases
- Fissure sealing

Ecosite Elements, Shades: White and Brown

Customization of composite restorations

Limitations of use

Application of the material is contraindicated if dry isolation or the recommended application technique is not possible.

Indication

- Loss of dental enamel due to caries, trauma, tooth wear or developmental disorders
- Primary and secondary prevention of caries

Contraindication

- Do not use directly on the exposed pulp.
- Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended user

Dentist

Notes for use

- For optimum polymerization, position the light-curing unit as close to the restorative material as possible and observe the recommendations relating to the light-curing time and the light intensity (see table). Observe the manufacturer's instructions for the light-curing unit used.
- Follow the manufacturer's instructions for other materials that are used with Ecosite Elements.
- If the material is applied to the patient's mouth with a single use applicator, the single use applicator must be used for this one patient only for hygienic reasons.

Recommended use

- Clean the tooth. Remove all cleaning residue with water spray. Determine the tooth shade using the Ecosite Shade Guide or VITA tooth shade. Select the suitable Ecosite Elements shade or shade combinations. Either choose a monochromatic Pure shade or combine it with a layer shade according to the overview. The instructions imply the use of a 0.5 mm-thick Ecosite Elements Layer.

Note: DMG recommends the use of a rubber dam.

- Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Ensure that the areas close to the pulp are protected.
- Condition the dental enamel and apply the bonding agent according to the manufacturer's instructions for the material used. DMG advises using the universal bonding agent Ecosite Bond.
- If necessary, apply a base of low-viscosity composite or compomer. We recommend using a matrix.
- Insert Ecosite Elements with a maximum layer thickness of 2 mm into the cavity. Apply taller restorations one layer at a time. Light-cure each layer separately.

Light intensity	≥ 1000 mW/cm ²	≥ 600 mW/cm ²
Exposure time	20 s	40 s
Max. layer thickness	2 mm	2 mm
Wavelength	420 – 480 nm	420 – 480 nm

Note: Additional buccal or lingual light-curing is advisable. Remove metal matrices beforehand.

- Removal of any excess material, contouring, and finishing are done using diamond finishing tools and flexible polishing disks of decreasing grain size. Ecosite Elements can be given a high shine using polishing paste together with polishing cups or disks.

Customization of composite restorations with W (White) and B (Brown)

- Select the color required for the desired color shade.
- Pour a sufficient amount of the selected shade onto a mixing palette.
- Apply the material to the composite with a fine brush in a thickness of ≤ 0.1 mm. Light-cure each layer separately.

Note: Always cover the customization layer with a composite layer.

Note: Mix the customization shades with one another (to obtain new tones), but never with other products.

Light intensity	≥ 1000 mW/cm ²	≥ 600 mW/cm ²
Exposure time	20 s	40 s
Max. layer thickness	0.1 mm	0.1 mm
Wavelength	420 – 480 nm	420 – 480 nm

Handling of the Safetip Dispenser Pro

- Pull the quick-lock mechanism all the way back.
- Insert the Safetip into the Safetip Dispenser Pro, and then guide the quick-lock mechanism back into place until the Safetip is clearly engaged.
- Turn the Safetip dispensing nozzle in the desired direction, remove the cap from the Safetip and dispense the material by pressing evenly. Do not exert too much force when pressing.
- After dispensing, pull the quick-lock mechanism back again and remove the Safetip.
- The Safetip Dispenser Pro should be cleaned and disinfected after use.

Disinfection of the Safetip Dispenser Pro

The Safetip Dispenser Pro is reusable after disinfection, although it should be replaced if visible signs of wear appear.

Disinfection can be carried out using a commercially available wiping or immersion disinfectant. DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute).

Interactions

- Materials containing eugenol, moisture and oily air may inhibit polymerization at the contact area. Avoid these conditions.
- The use of oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

Residual risks/side effects

There are no known side effects to date. The residual risk of hypersensitivity to components of the material cannot be ruled out.

Warnings/precautions

- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Reuse: Standard surface disinfection as required. No specific maintenance required.
- Serious incidents involving this product must be reported to the manufacturer and to the responsible registration authorities.

Product performance characteristics

Ecosite Elements in the shades: Bleach 1, B1, A1, A2, A3, A3.5, Enamel Bleach, Enamel Light, Enamel Medium, Enamel Dark complies with DIN EN ISO 4049: 2019 (Type 1, Class 2, Group 1); polymerization depth: ≥ 2 mm by using a light-curing time of 40 seconds with a light intensity of ≥ 600 mW/cm² or by using a light-curing time of 20 seconds with a light intensity of ≥ 1000 mW/cm²

Ecosite Elements in shades A2 Opaque and Incisal conforms to DIN EN ISO 4049:2019 (Type 2, Class 2, Group 1); polymerization depth ≥ 2 mm with a light-curing time of 40 s and a light intensity of ≥ 600 mW/cm² or a light-curing time of 20 s with a light intensity of ≥ 1000 mW/cm²

Storage/disposal

- Store in a dry place at 2 to 25°C (36 to 77°F)!
- Disposal must comply with national regulations.

Composition

Ecosite Elements BL1, B1, A1, A2, A3, A3.5, Enamel Bleach, Enamel Light, Enamel Medium, Enamel Dark:

Dental Glass, EBPADMA, bis-GMA, NPG2PODA, UDMA, SiO₂, additives
Inorganic filling materials: approx. 65 vol% (0.02 – 0.7 μm)

Ecosite Elements A2O, INC, W and B:

Dental glass, YbF₃, EBPADMA, bis-GMA, TEDMA, UDMA, SiO₂, additives
Inorganic filling materials: approx. 38 vol% (0.02 – 0.7 μm)

 Contains: 2-ethylhexyl 4-(dimethylamino)benzoate

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