





## Leading the way to a new standard

With over 1 billion composite restorations and a century of competence in dental materials, GC is leading the way to new standards in adhesive dentistry. With its proprietary **Dual H-Technology**, an advanced optimisation of the bonding both to the tooth and the composite, G2-BOND Universal is an **all-round universal 2-bottle system** combining what clinicians would expect from existing self-etch (SE) and etch-and-rinse (ER) Gold Standards, and even more:



- Your choice of preferred etching mode
- Superior durability and bond strength
- Superior resistance to marginal staining
- Superior dispensing and low technique sensitivity
- Multiple indications <sup>1</sup>

## Evaluations by clinicians worldwide resulted in over 1,500 restorations performed with G2-BOND Universal.<sup>2</sup> The results indicate:

- Remarkable sealing effectiveness with virtually no post-op sensitivity<sup>3</sup>
- 97% satisfaction in SE and ER etching modes
- 98% satisfaction with product handling and dispensing
- 89% satisfaction with the ergonomic bottle design



"The flexibility to bond in self-etch or etch-and-rinse mode, to seal and repair existing restorations is now applied to a 2-bottle adhesive." 4

<sup>&</sup>lt;sup>1</sup> Indications: Direct Composites I Intra-oral Repair I Immediate Dentin Sealing (IDS) I Treatment of hypersensitive teeth I Treatment of exposed root surfaces I Luting of indirect restorations (with G-CEM LinkForce™) I Post luting & Core Build-up (with GRADIA™ CORE & G-Premio BOND DCA).


<sup>&</sup>lt;sup>2</sup> Data on File


 $<sup>^{\</sup>rm 3}\,$  Results confirmed by independent study from University of Siena, Italy.

<sup>&</sup>lt;sup>4</sup> Prof Van Meerbeek, KU Leuven, Belgium.

## Leading the way to...

## ...Superior durability & higher bond strength







Source: Dr Mine, Osaka University, Japan

Source: Prof Van Meerbeek, KU Leuven, Belgium

G2-BOND Universal creates a **physically stronger**, **more durable** adhesive layer as the **Dual-H Technology** ensures a **higher bond strength immediately after placement**, with lower water sorption which **decreases the risk of degradation** for superior durability.

### ...Unmatched marginal quality






Source: Prof Tagami, Tokyo Medical and Dental University, Japan


Source: Prof Ferrari, University of Siena, Italy

A strong bond layer that assists the aesthetic transition from the composite to the tooth, **avoids gap formation and microleakage** to maintain marginal quality over the time.

### ...Shock absorbing layer



Source: (SEM picture) Prof Miyazaki, Nihon University, Japan (Film thickness) GC R&D, Japan, 2020, Data on File




Source: GC R&D, Japan, 2020, Data on File

The **optimal thickness** and **flexural strength** of G2-BOND Universal may act as a shock absorbing layer to **avoid debonding**. It functions **as a base** in large cavities and can overcome the polymerisation shrinkage, **compensating for the stress caused by the over layered composite.** Now clinicians can enjoy the benefits of a thick bonding layer in all etching modes.

## Leading the way to...

## ...Superior clinical handling

### **IDS** technique



Etching the tooth structure



Application of 1-PRIMER



Application of 2-BOND




After photopolymerisation, the sealed dentin is ready for digital impression.

#### **Class II restoration**



Selective enamel etching



Application of 1-PRIMER




Application of 2-BOND




After photopolymerisation, the restorative material is applied.

## ...a new 2-bottle bonding standard in all etching modes



G2-BOND Universal Clearfil SE Bond 2



G2-BOND Universal OptiBond FL

The spider graphs above show the all-round performance of G2-BOND Universal that leads the way to:

• Long-lasting restoration

of your patients!

- Virtually no post-operative sensitivity and aesthetic bonding layer
- Superior handling and dispensing

# For an all-round satisfaction and comfort



## G2-BOND Universal Dual H-Technology

Challenged to **improve on existing 2-bottle benchmarks**, GC researchers focused on strengthening the adhesive layer and significantly reducing the potential for water degradation over time. The new GC proprietary **Dual-H Technology provides a smooth transition from hydrophilic to hydrophobic characteristics** in the primed surface, which is then **further strengthened by a truly hydrophobic protective bond layer.** <sup>1</sup>

Hydrophilic primer wets and self-etches the tooth surface effectively.

Without hydrophilicity?

- Insufficient resin infiltration to dentin surface
- Higher risk for gap formation.



Highly hydrophobic bond layer results in strong bonding to composite.

Without hydrophobicity?

- Hydrolytic degradation of bond layer
- Increased risk for debonding.

#### All-round performance of a unique composition

- 4-MET for a reliable adhesion to enamel & dentin
- MDP for an excellent adhesion to enamel & dentin, zirconia, alumina and non-precious metals
- MDTP to bond to precious metals
- HEMA-free bonding system for a clinically proven outstanding durability <sup>2</sup>

G2-BOND Universal (0% HEMA)

Clearfil SE BOND 2 (10-30% HEMA)

Optibond FL (20-40% HEMA)

### Modern, intelligent and ergonomic design

G2-BOND Universal is instantly recognisable by its unique bottle covers

#### 1-PRIMER

- Bottle cover designed for low viscosity & volatile content to ensure controlled dispensing
- Prevent dripping and eliminate waste of liquid through ingenious cover design
- 300 applications per bottle
- Reusability with any new bottle

#### 2-BOND

- Designed for high viscosity liquids to reduce dispensing force
- Prevent dripping and eliminate waste of liquid through ingenious cover design
- 300 applications per bottle
- Reusability with any new bottle





Reduced contact point between cover and bottle





Rib designed to effectively press the bottle

<sup>&</sup>lt;sup>1</sup> GC R&D Date, Japan, 2020, Data on File.

<sup>&</sup>lt;sup>2</sup> Van Dijken J, et al, Dent Mater. (2013).

Order no.	Description	Content
10006819	G2-BOND Universal Bottle Kit	1x 1-PRIMER 5mL 1x 2-BOND 5mL accessoires
10006821	G2-BOND Universal 1-PRIMER Refill 5ml	1x Refill 5mL
10006822	G2-BOND Universal 2-BOND Refill	1 x Refill bottle 5mL



### **Related products**



G-ænial A'CHORD™ Advanced universal composite



**everX Flow™**Short-fibre reinforced flowable composite

G-ænial A'CHORD and everX Flow are trademarks of GC. Clearfil SE Bond 2 and OptiBond FL are not trademarks of GC.

#### GC EUROPE N.V.

http://europe.gc.dental

Head Office Researchpark, Haasrode-Leuven 1240 Interleuvenlaan 33, B-3001 Leuven Tel. +32.16.74.10.00 Fax.+32.16.40.48.32 info.gce@gc.dental

#### GC UNITED KINGDOM Ltd.

Coopers Court Newport Pagnell Buckinghamshire MK16 8JS United Kingdom Tel. +44.1908.218.999 Fax.+44.1908.218.900 info.uk@gc.dental http://europe.gc.dental/en-GB

