Icon Proximal

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English

Product description

Icon Proximal is an innovative product for the micro-invasive treatment of lesions in the proximal area.

Icon Proximal has been specially developed for the gentle treatment of incipient proximal caries.

Intended purpose (Fig. see “Lesion depth classification”)

Micro-invasive treatment of non-cavitated cariogenic lesions (lesion depth up to D1*)

*Radiographic lesion depth classification according to bitewing images.

Limitations of use

Do not use the material unless it is possible to ensure an adequate dry isolation and to protect the gingiva.

Indication

Non-cavitated cariogenic lesions (lesion depth up to D1*)

* Radiographic lesion depth classification according to bitewing images.

Contraindication

- Do not apply the material for deeper seated lesions (D2 – D3) (Fig. see “Lesion depth classification”) with cavitated enamel (enamel defects).
- Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended user

Dentist

Notes for use

- Light-curing units should emit at 450 nm and should be checked regularly.
- The light intensity should be a minimum of 800 mW/cm².
- Place the light unit as close to the material as possible.
- The etched surface must not be touched or contaminated with saliva or moisture during the treatment. If contamination occurs after drying, re-etch for approx. 10 seconds and repeat drying with Icon Prime.
- Do not use rubber dams made from thermoplastic elastomers such as Flexi-Dam (Coltène Whaledent/Hygenic).
- Additional proximal surfaces can be treated consecutively during one visit.
- The material included in the treatment unit contains sufficient material for treating two proximal lesions of one patient.
- Discard syringes after use.
- Do not press the plunger! Dispense material exclusively by turning the spindle.

Recommended use (Fig. see “Handling”)

1. Preparation

Before the start of treatment, clean the affected tooth and adjacent teeth. Remove all cleaning residues with water spray.

2. Isolation

Suitable measures for a complete isolation of the working area must be taken to make sure that the area is sufficiently dry.

Icon Etch can result in chemical burns to the oral mucosa if not used as directed. It is therefore essential to protect the gingiva completely by using a rubber dam or DMG MiniDam to prevent tissue damage.

3. Separation

- Introduce one of the enclosed dental wedges into the interproximal space (Fig. 1).
- For better access to the proximal region the wedge handle can be bent or detached with a twisting motion.
- In order to obtain sufficient separation of the teeth, the wedges must be seated at an adequate depth in the interproximal space.
- A separation of approx. 50 µm is recommended to ensure a successful treatment. This can be performed gradually: Insert the wedge to the point of resistance. Wait 3 to 5 seconds and gradually move the wedge deeper until sufficient separation is achieved.
- Leave the wedge in the proximal space during the entire treatment procedure.
- For tighter proximal areas, the use of an orthodontic rubber ring or other separation method prior to treatment will ensure adequate separation.

4. Application of Icon Etch

For effective treatment results the hypermineralization must be removed.

1. Screw the Proximal Tip onto the Icon Etch syringe, and introduce the application tip into the proximal area. Make sure that the green side of the Proximal Tip faces in the direction of the surface to be treated (Fig. 2). The material only dispenses on the green side of the Proximal Tip.
2. Apply a sufficient amount with a slight excess of Icon Etch to the lesion site by gently turning the spindle and then allow Icon Etch to act for 2 minutes, activating the material by occasionally moving the tip. About 1.5 to 2 turns of the syringe shaft correspond approximately to the amount to be applied. Remove excess material (Fig. 3).
3. Remove the application tip from the interproximal space. Suction off Icon Etch and rinse with water for at least 30 seconds (Fig. 4). Then dry using oil-free and water-free air.

5. Application of Icon Prime

For effective treatment results it is necessary to dry the lesion to its full depth.

1. Screw the application tip onto the Icon Prime syringe, turn the spindle to apply an ample amount of material onto the lesion and allow to take effect for 30 seconds (Fig. 5). Then dry carefully using oil-free and water-free air.

6. Application of Icon Infiltrant

Do not apply Icon Infiltrant to the tooth under direct operating light, as this may cause the material to cure prematurely.

1. Screw a new Proximal Tip onto the Icon Infiltrant syringe, and introduce the application tip into the proximal area. Make sure that the green side of the Proximal Tip faces in the direction of the surface to be treated. The material only dispenses on the green side of the Proximal Tip.
2. Apply an ample amount of Icon Infiltrant to the lesion site by turning the shaft. About 1.5 to 2 turns corresponds roughly to the amount to be applied.
3. Allow Icon Infiltrant to infiltrate for 3 minutes, occasionally activate the infiltration by moving the tip and top up if necessary (Fig. 6).
4. Remove the application tip from the interproximal space. Remove excess material with dental floss.
5. Light-cure Icon Infiltrant for a minimum of 40 seconds from all sides (Fig. 7).
6. Apply Icon Infiltrant again, leaving it to infiltrate for 1 minute, then remove any excess and light-cure it as described above.
7. Remove the wedge and rubber dam.
8. Use polishing strips for the surface finish.

Interactions

Do not use together with substances which inhibit polymerization. This especially applies to materials containing peroxide, such as bleaching gels.

Residual risks/side effects

- In some cases, contact allergies with similar composite products have been reported.
- Contact of Icon Etch with the oral mucosa causes a white surface coloration. This will subside after a few days.

Warnings/precautions

- Limit the application to the surface to be treated. Avoid contact of the etching gel with gingiva, exposed dentin, the tooth root and adjacent teeth.
- Use only the enclosed syringe and application tip to apply the etching gel in the patient’s mouth.

Caution! Icon Etch can cause chemical burns if not used properly.

- Avoid contact with soft oral tissue, eyes, and skin. In case of inadvertent contact rinse with copious amounts of water, and, if necessary, consult a physician.
- Discard the applicator tip after use. Replace syringe cap.
- Follow the manufacturer’s instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer (productsafety@mtcompanies.com) and to the responsible registration authorities.

Product performance characteristics

Penetration coefficient ≥ 100 cm/s

Polymerization depth ≥ 2 mm

Storage/disposal

- Store at 2 – 25 °C (36 – 77 °F) in a dry place.
- Protect from exposure to direct sunlight!
- Do not use after the expiration date!

Patient card + sticker (Fig. see “Patient card”)

This material is not radiopaque. For the documentation of the infiltrated surfaces and lesion depths, the enclosed patient card can be filled out at the time of the treatment and follow-up visits and be given to the patient.

- 1. Front side: Enter the name of the patient and add the office stamp.
- 2. Tooth diagram: Mark the treated area.
- 3. Record the treated tooth [1], the treated surface [2] (d = distal, m = mesial, v = vestibular (buccal, cervical, labial), l = lingual or palatal), lesion depth*, and the treatment date [3].
- 4. Follow-up visit [4]: Record the lesion depth* and date of the follow-up appointment. Annual check-ups are recommended.

*see “Lesion depth classification”

Composition

- Icon Etch: Water, hydrochloric acid, glycerol, SiO₂, additives
- Icon Prime: Ethanol, MDP
- Icon Infiltrant: TEDMA, TMPEOTA, additives

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Icon Proximal

Starter pack	
1 treatment unit, including: 1 Syringe @ 0.45 ml Icon Etch 1 Syringe @ 0.45 ml Icon Prime 1 Syringe @ 0.45 ml Icon Infiltrant 6 Icon Proximal Tips 1 Icon Capillary Tip 4 Icon Wedges	REF 220523
Pack	
2 treatment units, each including: 1 Syringe @ 0.45 ml Icon Etch 1 Syringe @ 0.45 ml Icon Prime 1 Syringe @ 0.45 ml Icon Infiltrant 6 Icon Proximal Tips 1 Icon Capillary Tip 4 Icon Wedges	REF 220524
Economy pack	
7 treatment units, each including: 1 Syringe @ 0.45 ml Icon Etch 1 Syringe @ 0.45 ml Icon Prime 1 Syringe @ 0.45 ml Icon Infiltrant 6 Icon Proximal Tips 1 Icon Capillary Tip 4 Icon Wedges	REF 220525

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