# Bluephase® Style

The curing light

# The smallest LED for every use







## **Bluephase**: The LED for every use ...



True innovation proves itself in use. A perfect example of this is the Bluephase curing light with the specially developed Polywave® LED, which is capable of polymerizing all current photoinitiators and materials. This puts the Bluephase light on an equal footing with the halogen lights, the previous industry standard.

Polywave

Bluephase sets new standards in every respect. It is not only the best-selling LED polymerization light in Europe, but it is also frequently used in clinical and scientific studies¹ and has been ranked top light-curing unit by independent and renowned test institutes.

Grounded in proven concepts and open to further development: These principles are embodied in the Bluephase Style. Whilst offering a performance comparable to that of the existing Bluephase light, the new **Bluephase Style** features a compact and ergonomic shape that **fits the hands of women and men** even more comfortably.



Bluephase: now even smaller

reddot design award

<sup>&</sup>lt;sup>1</sup> IADR Barcelona 2010 Mentions in abstracts

# ... and now comfortable to fit every hand



Conventional LED devices are available either as a pistol or a pen. Bluephase Style combines the two: It can be held like a **pen** or like a **pistol**. Given the light's low weight and balanced weight distribution, the strain exerted on the hand and arm is reduced, regardless of the size of the hand.

The light probe of Bluephase Style clearly stands out from all the other known light probes. It is **shortened** and can be rotated by 360 degrees, enabling easy access to all tooth surfaces **without extreme opening of the mouth**.





Two-button operation



Bluephase Style can be held like a pen or like a pistol.



Competitive light



Even difficult-to-reach areas can be easily accessed. (Dr Eduardo Mahn, Santiago, Chile)

#### **Every material** due to Polywave® LED

#### Every indication due to continuous cooling

#### Every time due to Click & Cure

1,100 mW/cm<sup>2</sup>


The capability to cure all dental materials is related to the light emitted. In contrast to the LED lights of the second generation, the **Polywave LED** curing lights achieve an optimal broadband spectrum of 385 to 515 nm, similar to the spectrum of halogen lights that served as their model. Therefore, Bluephase Style is suitable for polymerizing all current photoinitiators and materials without any restrictions.

Because Bluephase Style incorporates particularly energy-efficient LED technology, it produces less heat than an LED light providing the same output. As a result, the need for a **fan** has been **eliminated**. That means you can just get on with your work even when using the light in continuous operation to place extensive indirect restorations.

Cordless operation ensures maximum mobility in the dental practice. The proven Click & Cure function enables you to avoid irritating waiting times in case the battery has run out. With just one click you can connect the handpiece to the power cord of the charging base.

**Easy operation** Two-button operation

#### Wavelength range



Source: R. Price, Dalhouse University Halifax, 2011

- \* These trade names are not registered trademarks of Ivoclar Vivadent AG.
- <sup>1</sup> LED light of the 2<sup>nd</sup> generation
- <sup>2</sup> LED light of the 3<sup>rd</sup> generation

"LED curing lights with polywave technology should be preferred over single peak devices." Price R, Fahey J, Felix C: Knoop hardness of five different composites cured with single peak and polywave LED curing lights. Quintessence Int. 2010; 41; e181-e191

#### Photoinitiator: camphorquinone

Photoinitiator: Lucirin TPO Bluephase Style with Polywave LED

Elipar S10 \*1

SmartLite PS \*1

Demi Plus \*1





This is how it works: Turn around charging base, remove power cord, attach to handpiece, resume work as usual.

**Inductive charging** enables contactless charging of the battery.

Charging base with battery charge indicator

The illuminated ring appears black if the battery is fully

charged.



### Bluephase® Style cures effectively even from a distance

Given the large **10-mm diameter** of the light probe, the light completely covers large cavities. Time-consuming multiple curing cycles for MOD restorations are thus no longer necessary.

Special optical technology enables the intensive light of Bluephase Style to penetrate deep into the light-curing material and ensures that an exceptionally high light intensity is available in critical areas in particular. A reliable cure can be achieved even when the material is light-cured from a long distance, for instance in proximal boxes.

Because of its high light intensity, Bluephase Style allows large curing depths to be achieved with curing times as short as 10 seconds for composites and adhesives.

#### Bluephase Style




Competitive light




The 10-mm light probe allows MOD fillings to be irradiated only once, saving valuable time.

## Transmission energy at a distance of 4 mm



Source: Excerpt R. Price, Dalhouse University Halifax, 2011 (study commissioned by Ivoclar Vivadent AG)

#### Curing depth of Tetric EvoCeram® A3



Curing time with Bluephase Style: 10 seconds; curing time with other devices: according to manufacturer's directions.

Source: R&D Ivoclar Vivadent Schaan, 2010



<sup>\*</sup> These trade names are not registered trademarks of Ivoclar Vivadent AG.

# **Bluephase** Style is available in three different colours



Bluephase Style is an eye catcher for your dental office. Available in a pink, blue and grey version, the small, handy polymerization light adds colour to daily office routines. And it fits perfectly into every practice environment.





# Bluephase® Style at a glance

# Technical data

Bluephase Style	1,100 mW/cm² ±10 %
Every hand (ergonomic design)	✓
Every material (wavelength range)	<b>✓</b> (385–515 nm)
Every indication (continuous operation for at least 10 min)	✓
Every time (Click & Cure: optional mains operation)	✓ (lithium-polymer battery)
Curing time for selected composite materials 2 mm   Tetric EvoCeram*/IPS Empress Direct* 4 mm   Tetric EvoCeram* Bulk Fill	10 s 10 s
Power supply	Contactless charging via lithium-polymer battery, capacity: approx. 20 min; charging time: approx. 2h
Weight of handpiece	120 g (incl. battery and light probe)
Dimensions of handpiece (without light probe)	L = 180 mm, W = 30 mm, H = 30 mm
Weight of charging base	195 g
Dimensions of charging base	D = 125 mm, H = 70 mm
Warranty	3 years (battery 1 year)

Every material, every indication, every time – and now for every hand. Only this unique combination gives you the licence to cure.

# Delivery forms/accessories

	Art. no.
Bluephase Style 100–240 V grey	635 153
Bluephase Style 100-240 V blue	642 513
Bluephase Style 100–240 V pink	642 514
Handpiece, battery, light probe 10 mm (black), charging base, power cord and power pack, anti-glare cones (3x), protective sleeves (50x), Operating Instructions	
Light probe 10 mm, black	636 240
Light probe Pin-Point 6 > 2 mm, black	636 241
Protective sleeves, 50x	636 239
Anti-glare cone	551 756
Anti-glare shield	592 496
Bluephase Style handpiece grey	637 916
Bluephase Style handpiece blue	643 296
Bluephase Style handpiece pink	643 307
Handpiece, battery, light probe 10 mm (black )	



#### **Direct Restoratives**



Bluephase® Style forms a part of the "Direct Restoratives" product category. The products of this category cover the procedure involved in the direct restoration of teeth – from preparation to restoration care. The products are optimally coordinated with each other and enable successful processing and application.



#### THESE ARE FURTHER PRODUCTS OF THIS CATEGORY:

#### Tetric EvoCeram® Bulk Fill

The bulk-fill composite



#### Posterior efficiency, zero compromise

- Fill with up to 4 mm bulk increments due to lvocerin®, the patented light initiator
- Sculpt and contour with ease as a result of the material's smooth consistency
- Then light-cure ... and the esthetic restoration is done!

#### **OptraGate**°

Latex-free lip and cheek retractor



#### The gentle solution for a better view

- Considerably enlarged access to the treatment field
- Enhanced patient comfort due to three-dimensional flexibility
- Universal range of indications

Would you like to know more about the products of the "Direct Restoratives" category?

Simply get in touch with your contact person at Ivoclar Vivadent or visit www.ivoclarvivadent.com for more information.

Ivoclar Vivadent AG

Rendererstr. 2 FL-9494 Schaan Liechtenstein Tel. +423 / 235 35 35 Fax +423 / 235 33 60

www.ivoclarvivadent.com

