Oral and Microsurgical Cutting Instrument

Dental Surgical Scissors, Suture Scissors, Periodontal Knife , Scalpel Handle

en-English

Attention

The operator/user must carefully read and understand this manual thoroughly to keep the product performance durable and reliable for defined life expectancy. After opening the packaging, first of all it is necessary to check the component(s) against the standard configuration. Check that these are all present and in perfect condition.

Intended Use:

Oral and microsurgical cutting instruments are intended to be used for precise cutting and incision of soft tissues and other materials during dental and oral surgical procedures. They assist in achieving controlled tissue manipulation, contouring, and trimming.

The Instrument Specific Intended Use is given below:

Dental Surgical Scissors:

Intended to be used for cutting soft tissue and materials during oral and periodontal surgical procedures, ensuring precise incisions and controlled tissue trimming.

Suture Scissors:

Intended to be used for cutting and removing sutures during post-surgical care, ensuring clean cuts without fraying or damaging surrounding tissue(s).

Periodontal Knife:

Intended to be used for making precise incisions in gingival and periodontal tissues during periodontal and oral surgical procedures, assisting in contouring and reshaping tissues as part of periodontal treatment.

Scalpel Handle:

Intended to be used for Holding and stabilizing interchangeable scalpel blades during dental and oral surgical procedures, allowing for precision and control when making incisions in soft tissue.

Our Oral and microsurgical cutting instruments are intended to be used by qualified dental/healthcare professionals only in a professional dental/healthcare facility.

Patient population:

These instruments are intended for use on adult and child patients undergoing oral and microsurgical procedures, including soft tissue cutting, suturing, and periodontal interventions. The appropriate instrument type and size must be selected based on the patient's oral anatomy, surgical site, and procedural requirements.

For Available Configurations/Models please visit the Oral & Microsurgery section in our Catalogue

Features

Dental Surgical Scissors & Suture Scissors:

- Sharp, fine-edged blades ensure precise cutting of soft tissues and sutures without fraying or damage.
- Available in curved and straight blade designs for versatility in different surgical applications.
- Ergonomic handles provide a secure grip and reduce hand fatigue, ensuring precise control in confined spaces.

Periodontal Knife:

- Sharp, fine-edged blades ensure precise cutting of soft tissues and sutures without fraying or damage.
- Available in curved and straight blade designs for versatility in different surgical applications.
- Ergonomic handles provide a secure grip and reduce hand fatigue, ensuring precise control in confined spaces.

Scalpel Handle:

- Designed to securely hold various blade sizes for specific procedures.

- Available in straight or angled designs for improved access to confined areas.
- Manufactured from high-quality stainless steel, ensuring durability, corrosion resistance, and full sterilization compatibility.

Operating instructions

 These instruments are non-sterile and must be cleaned and sterilized before first use.

 Always check if the reprocessed instruments are in perfect condition and do not use the instruments if any kind of deterioration, damage or deformation is observed.

 Use a firm but controlled grip to maintain precision and stability. Avoid excessive force to prevent unnecessary trauma to soft tissues and surgical sites. Ensure all cutting instruments are sharpened (as and if needed) before use to maintain efficiency and precision.

For Dental Surgical & Suture Scissors:

- Select the appropriate scissors based on the procedure.
- Position the blades carefully to ensure a clean and precise cut.
- Use controlled pressure to avoid excessive force that may cause tissue trauma.

 Always check the hinge movement before use. Stiff or loose hinges may affect performance and should be addressed before the procedure.

For Periodontal Knife:

- Ensure proper blade selection for the intended surgical procedure.
- Use gentle, controlled strokes to make precise incisions while minimizing tissue damage.
- Maintain a steady grip and correct angulation to ensure optimal tissue contouring.

For Scalpel Handle:

- Securely attach the appropriate scalpel blade to the handle before use.
- Ensure the blade is properly locked in place to prevent slippage during incisions.
- Hold the handle with a stable grip, applying controlled pressure for accurate cutting.

 When handling scalpel handles, ensure blades are securely attached and properly disposed off after use to prevent injury.

Adverse Event

Any serious incident that has occurred concerning the device should be reported to the manufacturer and to the respective competent authority of the Member State in which the incident has occurred.

Cleaning & Sterilization

Cleaning

 Clean the device with care after each application to avoid the drying of contaminants.

Initial Cleaning

1. Remove any visible debris from the instrument using moisten a cotton swab or cloth with isopropyl or ethyl alcohol and gently wiping the surface of the product.
2. Use a soft brush and enzymatic cleaner to remove contaminants thoroughly.

Manual Cleaning

3. Prepare a solution of neutral pH enzymatic detergent according to the manufacturer's instructions.
4. Immerse in a pre-soak enzymatic cleaner solution
5. Thoroughly clean all surfaces of the Instrument using a soft brush or cloth. Pay close attention to any crevices or hard-to-reach areas.
6. Rinse thoroughly under lukewarm running tap water for a minimum of 30 seconds to remove all traces of detergent.

Automated Cleaning

1. Always clean in an appropriate box or cassette and ensure that they are properly positioned within the unit.
2. Always follow the equipment manufacturer's instructions and ensure the instruments are compatible with the cleaning system.
3. Visually inspect for cleanliness, If any contamination is visible, repeat the cleaning steps until required cleanliness is achieved.
4. Dry with clean, lint free cloth or filtered compressed air until there is no visible moisture. Follow STEAM STERILIZATION PROCEDURE.

⚠ For scissors, ensure that the hinge mechanism is fully cleaned to prevent residue buildup that could affect movement.

⚠ Ultrasonic cleaning may be performed in accordance with hospital or facility protocols. Ensure proper rinsing and drying after the process to prevent residue buildup.

Sterilization

Note: Ensure that the instruments are completely dry before sterilization.

Steam Autoclaving

1. Place the instrument in an approved sterilization pouch or wrap.
2. Sterilize using a steam autoclave following a validated sterilization cycle according to ISO 17665-1 and your facility's protocols. A typical recommended cycle is 132-134°C (270-273°F) for 4 minutes.
3. Verify no residual moisture remains after sterilization to prevent corrosion and ensure longevity.
4. Ensure that the hinges are lubricated with an appropriate medical-grade lubricant after sterilization to maintain smooth operation.
5. Allow the instruments to cool and dry completely before handling. Store in a clean, dry, and controlled environment until ready for use.
6. Follow the pouch/wrap manufacturer instructions for storage conditions and maximum storage time.
7. For scissors, periodically check the hinge movement to ensure continued smooth function.

Important Notes

- Always Inspect the instruments for any signs of deterioration, breakage and bent etc. after each sterilization cycle. Discard if damage is observed.
- Always follow your facility's infection control protocols.
- Refer to your autoclave manufacturer's instructions for proper operation and maintenance. The Sterilization equipment must be validated by the hospital and or sterilization equipment manufacturers.
- The above-mentioned cleaning and sterilization guidelines, provided by manufacturer are intended as procedures compatible with specific materials. Sterilization must be performed according to the Hospital/Clinic approved protocol.

- The responsibility for ensuring proper sterilization of instruments lies with the user. The effectiveness of sterilization depends on validated cleaning, packaging, and sterilization procedures carried out at the facility.

Disposal

Our Oral and Microsurgical Cutting Instrument: Dental Surgical Scissors, Suture Scissors, Periodontal Knife and Scalpel Handle contain no hazardous materials. However, these must be cleaned/sterilized as per defined procedures before disposition. Please follow the local and national regulations or healthcare facilities' defined disposal and waste management policies.

Disclaimer

The products must be used, reprocessed, and maintained strictly in accordance with the instructions provided above. Any deviation from these guidelines by the dental professional or user is undertaken at the user's sole risk. The Manufacturer will not accept any requests for refunds or exchanges under warranty for products that have not been handled and reprocessed in compliance with these instructions.

Explanation of utilized symbols

	Medical Device Compliant with Regulation (EU) 2017/745
	Consult instructions for use
	Medical Device
	Product Code/Catalogue Number
	Model/Variant
	Lot Number
	The Device is supplied non-sterile.
	Sterilizable in a steam sterilizer (autoclave) at 132-134°C (270-273°F)
	Manufacturing Date
	Manufacturer
	Authorized representative in the European Community

	Importer
	Keep in a cool, dry place
	Caution, read instructions

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