

TePe PlaqSearch™

Disclosing tablet	Ingredients: Hydrated Dextrates, Magnesium Stearate, Sodium Starch Glycolate, Aroma, CI 42090, CI 45410.
Pack	10 tablets in a plastic strip/consumer pack. 12 consumer packs/display carton. 250 tablets in a multi-pack (10 tablets/plastic strip, 25 plastic strips/carton).
Recycling	Plastic strip in container for plastic, box in container for cardboard. The carton should be placed in container for corrugated cardboard.
Product features	 A disclosing product for quick and simple identification of plaque. Plaque is difficult to see, especially between the teeth. Differentiates new (red) and old (blue) plaque. This is important as the older plaque can be more harmful. Improves toothbrushing. Locates the areas where plaque has accumulated and where brushing needs to be directed. Avoid contact with clothes and furnishings.

- Children under 12 years should be supervised by an adult.

