




Urban General Surgery Kit



Instruments are placed in a surgical cassette according to the order in which they will be used during the surgery. These instruments are sterilized in the surgical cassette before the surgery. The instruments are used out of the cassette and are kept in order during the entire surgical procedure.

Urban General Surgery Kit

| | | | |
|--|-------------------|---|-------------------|
| IMS™ Infinity Series™ DIN Surgical Cassette for 18 Instruments with an accessory area + 4 clips (2 x I1000) for hinged instruments | IMEDIN18_S | Scissors Goldman-Fox #16 curved serrated, SuperCut 12,5cm | S16SC |
| Mirror Handle hdl #6 US-Thread | MH6 | Tissue Pliers De Bakey PermaSharp™ straight 15 cm | TP5090 |
| High Definition CS Mirror #5 single-sided, Cone socket, 1pc/pack | MIR5HD | Cheek Retractor Minnesota | CRM |
| Tweezers #17 | DP17 | Friescio Mouth prop #54, adult | MP54 |
| Probe UNC #15/12 hdl #6 | PCNC15/126 | Hemostat Halsted-Mosquito #4 straight, 12 cm | H4 |
| Scalpel Handle #5E, European handle | 10-130-05E | IMS™ 1/2 Bur Cushion Short Lid (6 Burs) | IMS-1372SH |
| Mini Me Microsurgical Periosteal, Black Line | PFIWDS1MKX | Periotome #1 Black Line, posterior | PT1X |
| Curette Columbia #4R/4L, #9 Handle, Everedge 2.0™ | SC4R/4L9E2 | Urban Chisel | CURBAN1 |
| Periosteal Prichard #3, Black Line | PPR3X | | |
| Scaler Ball #B-2/3 hdl #6 | SB2/36 | | |
| Tissue Pliers Micro Cooley hdl #SinusLine straight modified 18cm | TPSLCOSM | | |
| Needle Holder Simion #5024R, round, PermaSharp™, straight, 18cm | NH5024SIM | | |

CHOOSE THE COLOR OF CASSETTES.

Just add the color coding number to the last part of the partcode (replace " " with the number).

| | | | | | |
|--------|---|----------|------------|---|-----------|
| Red | ● | 0 | Lavender | ● | 6 |
| Grey | ● | 1 | Black | ● | 7 |
| White | ○ | 2 | Blue | ● | 8 |
| Orange | ● | 3 | Green | ● | 9 |
| Purple | ● | 4 | Ocean Blue | ● | 11 |
| Yellow | ● | 5 | | | |



Learn more about our Key Opinion Leader **Istvan Urban**.