

Cordless Obturation System

Fast-Pack





Fast-Fill



Fast-Fill

- Useable cartridge tip, economical design
- GP pellet refill into device
Easy maintenance



- Motorized gutta-percha
Push & stop, better flow controlling

Push Speed
Mid

- 3 ranges pushing speed, low, mid, high



- Powerful battery


T1 150°C
GP 100%

Temperature
160°C


- Digital displaying of rest gutta percha
- Real-time temperature displaying
Set from 100~200°C




Fast-Pack



- Ring switch, 360° touchable



- 4-6mm tip heating , safe and fast



- Real-time temperature displaying

- Customized temperature from 90~250°C



www.eighteeth.com

