

VDW.ULTRA®

High performance endo ultrasonic device

- Endo focus with clearly marked endo power ranges in large scale
- Correct power setting within the power range allows optimal balance between efficient and gentle operation
- Patented piezo technology



VDW.ULTRA® _____

- Incl. ultrasonic handpiece
- Incl. 1 ultrasonic endo tip

V04 1400 000 501

Ultrasonic Handpiece _____

V04 1401 000 000



VDW Endo Ultrasonic Tips

CAVI tips (diamond-coated)

CAVI 1-BD (ball)

- Access refinement

CAVI 2-D (swan neck, optimised visibility)

CAVI 3-D (standard)

Blister pack of 1 tip

CAVI 1-BD	V04 1410 000 001
CAVI 2-D	V04 1410 000 002
CAVI 3-D	V04 1410 000 003



CAVI 1



CAVI 2



CAVI 3

REDO tips

REDO 1 (short)

REDO 2 (medium, special Ti-Nb alloy)

REDO 5 (long)

- Retreatments

Blister pack of 1 tip

REDO 1	V04 1411 000 001
REDO 2	V04 1411 000 002
REDO 5	V04 1411 000 005



REDO 1



REDO 2



REDO 5

MAXI tip

MAXI MPR

- Dislodging metal posts at high power level

Blister pack of 1 tip

MAXI MPR	V04 1412 000 000
----------	------------------



MAXI

IRRI files

IRRI S (smooth wire)

- Passive ultrasonic irrigation (Pull)

IRRI K (K-type)

- Exposure of delicate anatomical structures

Blister pack of 4 files Length/ISO

IRRI S	21/25	V04 1460 021 025
IRRI S	25/25	V04 1460 025 025
IRRI K	21/15	V04 1463 021 015
IRRI K	21/25	V04 1463 021 025
IRRI K	25/15	V04 1463 025 015
IRRI K	25/25	V04 1463 025 025



IRRI S



IRRI K

VDW.ULTRA® is a distribution good.

Legal manufactory is SATELEC, a Company of ACTEON Group - 17, avenue Gustave Eiffel - BP 30216 - 33708 MERIGNAC Cedex - FRANCE