


PRODUCT	Hawe Matrices Super Thin																					
REFERENCE	<table border="1"> <thead> <tr> <th colspan="3">Hawe Matrices Super Thin</th> </tr> </thead> <tbody> <tr> <td rowspan="8">Refill pack of 30</td> <td>Hawe Matrices Super Thin</td> <td>1201</td> </tr> <tr> <td>Hawe Matrices Super Thin</td> <td>1202</td> </tr> <tr> <td>Hawe Matrices Super Thin</td> <td>1203</td> </tr> <tr> <td>Hawe Matrices Super Thin</td> <td>1204</td> </tr> <tr> <td>Hawe Matrices Super Thin</td> <td>1205</td> </tr> <tr> <td>Hawe Matrices Super Thin</td> <td>1206</td> </tr> <tr> <td>Hawe Matrices Super Thin</td> <td>1207</td> </tr> <tr> <td>Hawe Matrices Super Thin</td> <td>1208</td> </tr> </tbody> </table>			Hawe Matrices Super Thin			Refill pack of 30	Hawe Matrices Super Thin	1201	Hawe Matrices Super Thin	1202	Hawe Matrices Super Thin	1203	Hawe Matrices Super Thin	1204	Hawe Matrices Super Thin	1205	Hawe Matrices Super Thin	1206	Hawe Matrices Super Thin	1207	Hawe Matrices Super Thin
Hawe Matrices Super Thin																						
Refill pack of 30	Hawe Matrices Super Thin	1201																				
	Hawe Matrices Super Thin	1202																				
	Hawe Matrices Super Thin	1203																				
	Hawe Matrices Super Thin	1204																				
	Hawe Matrices Super Thin	1205																				
	Hawe Matrices Super Thin	1206																				
	Hawe Matrices Super Thin	1207																				
	Hawe Matrices Super Thin	1208																				
COMPANY/ Manufacturer	KerrHawe SA																					
CLASSIFICATION	CLASS I MEDICAL DEVICE																					
PRODUCT DESCRIPTION	Stainless steel matrices																					
SPECIFICATIONS	<p>1201:</p> <ul style="list-style-type: none"> • Material: Stainless steel 1.4310 • Overall dimensions: 50.5 mm x 10.5 mm • Thickness: 0.038 mm <p>1202:</p> <ul style="list-style-type: none"> • Material: Stainless steel 1.4310 • Overall dimensions: 50.2 mm x 10.4 mm • Thickness: 0.038 mm <p>1203:</p> <ul style="list-style-type: none"> • Material: Stainless steel 1.4310 • Overall dimensions: 50.7 mm x 11 mm • Thickness: 0.038 mm <p>1204:</p> <ul style="list-style-type: none"> • Material: Stainless steel 1.4310 • Overall dimensions: 51.8 mm x 11.1 mm • Thickness: 0.038 mm <p>1205:</p> <ul style="list-style-type: none"> • Material: Stainless steel 1.4310 • Overall dimensions: 64.4 mm x 14.6 mm • Thickness: 0.038 mm <p>1206:</p>																					

	<ul style="list-style-type: none"> • Material: Stainless steel 1.4310 • Overall dimensions: 66.8 mm x 11.1 mm • Thickness: 0.038 mm <p>1207:</p> <ul style="list-style-type: none"> • Material: Stainless steel 1.4310 • Overall dimensions: 64.5 mm x 12.2 mm • Thickness: 0.038 mm <p>1208:</p> <ul style="list-style-type: none"> • Material: Stainless steel 1.4310 • Overall dimensions: 65.6 mm x 13.2 mm • Thickness: 0.038 mm
PACKAGING	Refill pack of 30 pcs: 5.8 cm width 8.6 cm length, 2.9 cm height
SHELF LIFE	1800 days
INTENDED OF USE	Use for both molar and premolar restorations
SECURITY	Only professional use
TRANSPORT AND STORAGE CONDITIONS	Ambient temperature