

# DIATECH® ECO LINE POLISHERS FOR COMPOSITES

The perfect moment needs preparation




## Advantages at a glance

- Reliable polishers for all types of composites
- Polishers are flexible and adapt to tooth shape
- Easy contouring and adjusting with the coarse Removers (white)
- Yellow Smoothers reduce surface roughness
- Grey Shiners leave the patient with an esthetic shine

ECO Line Polishers for Composite Kit  
Order No. 60023468

All shapes and grit sizes  
combined in one kit.



# ORDER INFORMATION

## REMOVERS



Shape	CUP	POINT	DISC
Polisher type	RE8005RA	RE8008RA	RE8009RA
Ø 1/10mm	065	050	100
Head L in mm	10	10	1
Order No. (Box of 12)	60023469	60023470	60023471

## SMOOTHERS




Shape	CUP	POINT	DISC
Polisher type	SM8005RA	SM8008RA	SM8009RA
Ø 1/10mm	065	050	100
Head L in mm	10	10	1
Order No. (Box of 12)	60023472	60023473	60023474

## SHINERS



Shape	CUP	POINT	DISC
Polisher type	SH8005RA	SH8008RA	SH8009RA
Ø 1/10mm	065	050	100
Head L in mm	10	10	1
Order No. (Box of 12)	60023475	60023476	60023477



Cup shapes are for polishing cusp slopes of posterior teeth, proximal line angle areas and buccal, lingual, and palatal surfaces of anterior restorations.



Points are ideal for polishing both the cervical area and the occlusal grooves and reliefs.



Disc shapes are used on incisors, large filling surfaces, incisal edges and interproximal areas.

© COLTENE – www.coltene.com

**UK** **Coltene/Whaledent Ltd.**  
The President Suite  
Kendal House, Victoria Way  
Burgess Hill, West Sussex  
RH15 9NF / UK  
T +44 1444 235486  
F +44 1444 870640  
info.uk@coltene.com

**CH** **Coltene/Whaledent AG**  
Feldwiesenstrasse 20  
9450 Altstätten / Switzerland  
T +41 71 757 53 00  
F +41 71 757 53 01  
info.ch@coltene.com



**COLTENE**