

Product information

SENDOLINE AB

K-Reamer

Compared to the K-File the K-Reamer has an open spiral design with three cutting edges, a 60° tip and a 120° cutting angle.

It is designed to effectively enlarge the root canals. The instrument has an excellent debris transportation.

The K-Reamer should be used with a screwing movement.

The K-Reamer is also available in nickel titanium name [NiTi-K Reamer](#).

ISO/Sizes: 06-140

Length: 21 mm, 25 mm, 28 mm, 31 mm

Package: Individual or assorted (in sizes 15-40 or 45-80) in blister. All packages come with 6 pcs/blister

Material: Operating part: Stainless steel. Handle: Plastic. Stopper: Silicone.

