

# OptraGate<sup>®</sup>

The latex-free lip and cheek retractor



The gentle solution  
for a better view

  
**ivoclar**  
**vivadent**<sup>®</sup>  
passion vision innovation

# OptraGate®

The latex-free lip and cheek retractor

OptraGate® retracts the lips and cheeks evenly and gently. Therefore, the need for using a dental mirror to assist in retracting soft tissue is eliminated. The focus can be placed on the actual work at hand. This heightens the efficiency and enhances the quality of the treatment.

**Easy access to the cavity:**

- Even retraction of the lips and cheeks
- Easily accessible, manageable treatment field

**Stress-free treatment:**

- Easy placement
- Relieves the dental team and the dental practitioner



Wide range of indications



Diagnostics



Intraoral scans



Prosthetics



Tooth whitening



Professional care



Restorative treatment

# OptraGate®

## Simplified isolation

In addition to the success of the dental treatment, patient satisfaction is an important consideration for clinicians. Unfortunately, it is often impossible to achieve a completely dry treatment field and ensure a full view of the oral cavity at the same time. In combination with other auxiliary aids such as cotton rolls and saliva ejectors, OptraGate provides a solution in many cases, as it allows adequate isolation and excellent visibility to be achieved whilst offering high wearing comfort.

### High patient comfort:

- Three-dimensionally flexible for a very good fit
- Soft material with recesses to accommodate the lip frenula

### Moisture control made easy:

- Quick and effective isolation
- More time to focus on the actual treatment procedure



Professional care in children



Restorative treatment in children



Pediatric dentistry



Orthodontics



Diagnostics in children



# OptraGate®

In blue and in pink

## Now in colour!

Due to its soft and flexible quality, OptraGate also fits children very comfortably. OptraGate is available in two new colours, which are designed to please children and increase their compliance during treatment procedures. As a result, the tiresome retraction of soft tissue is turned into an exciting experience for adults and children alike.

### Good tolerance:

- Latex-free
- Gentle to the oral tissues
- Flexible fit

### Making treatment fun:

- Two colours
- Enhanced patient compliance



OptraGate® – Easy placement, easy procedure. See for yourself!

[www.ivoclarvivadent.com/OptragateFlowchart](http://www.ivoclarvivadent.com/OptragateFlowchart)



## Easy placement



A relaxed mouth makes placement easier



Insertion in the right corner of the mouth



Insertion in the left corner of the mouth



Placement behind the lower lip



Placement behind the upper lip



Optimally placed OptraGate

# Delivery forms

## OptraGate white



**OptraGate Assortment** 577275WW  
40 x size Regular  
40 x size Small

**OptraGate Regular Refill** 590850WW  
80 x size Regular

**OptraGate Small Refill** 590851WW  
80 x size Small

**OptraGate Junior Refill** 591451WW  
80 x size Junior

**NEW**

**OptraGate Regular Refill** 683377WW  
40 x size Regular

**NEW**

**OptraGate Small Refill** 683378WW  
40 x size Small

**NEW**


**OptraGate Junior Refill** 683376WW  
40 x size Junior

## OptraGate blue and pink

**OptraGate Starter Kit/blue and pink** 669087WW  
4 x size Junior blue  
4 x size Junior pink  
4 x size Small blue  
4 x size Small pink

**OptraGate Small Assortment/blue and pink** 669088WW  
20 x size Small blue  
20 x size Small pink

**OptraGate Junior Assortment/blue and pink** 669089WW  
20 x size Junior blue  
20 x size Junior pink





## Direct Restoratives

OptraGate® forms a part of the "Direct Restoratives" product category. The products of this category cover the procedure involved in the direct restoration of teeth – from preparation to restoration care. The products are optimally coordinated with each other and enable successful processing and application.



### THESE ARE FURTHER PRODUCTS OF THIS CATEGORY:

#### Tetric EvoCeram® Bulk Fill

The future of composite technology.  
Now available.



The fast posterior composite

- Bulk placement up to 4 mm due to Ivocerin®, the new light initiator
- Low shrinkage and low shrinkage stress for superior margins
- Sculptable consistency, extended working time under operatory light

#### Fluor Protector S

The protective fluoride varnish



Homogeneous protective varnish  
in a dispensing tube

- Immediate fluoridation of enamel
- High-yielding fluoride depot
- Mild taste

Would you like to know more about the products of the "Direct Restoratives" category? Simply get in touch with your contact person at Ivoclar Vivadent or visit [www.ivoclarvivadent.com](http://www.ivoclarvivadent.com) for more information.

Hersteller und Vertrieb  
**Ivoclar Vivadent AG**  
 Bendererstr. 2  
 9494 Schaan  
 Liechtenstein  
 Tel. +423 235 35 35  
 Fax +423 235 33 60  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

Darstellungen und Angaben enthalten keine  
Zusicherung von Eigenschaften.

Gedruckt in Liechtenstein  
© Ivoclar Vivadent AG, Schaan/Liechtenstein  
665085/2017-02-27/EN

**ivoclar**  
**vivadent**  
 passion vision innovation