

# Getting Started with the Versah® Guided Surgery System

1.

## Determine the Following:

- Tooth #
- Implant Diameter

### Versah® Guided Surgery System Price Breakdown

| Item                       | Unit Price |
|----------------------------|------------|
| G-Stop™ Vertical Gauge ... | \$24 each  |
| G-Stop™ Key .....          | \$69 each  |
| C-Guide™ Sleeve.....       | \$19 each  |
| G-Stop™ Holder .....       | \$95       |

2.

## Plan the Case:

- Take a CBCT Scan, Dicom File, and impressions of the site and the opposing jaw.
- Contact a Versah® Certified Lab to begin the process of printing a custom C-Guide™ Sleeve.

## Versah Certified Labs

A full list of certified labs available at [versah.com/guided-surgery-system/](http://versah.com/guided-surgery-system/)

### Dorsey Dental

[www.dorseydentallab.com](http://www.dorseydentallab.com)

### Cinzara

[www.cinzara.com](http://www.cinzara.com)

### ROE Dental Lab

[www.roedentallab.com](http://www.roedentallab.com)

### 3D Diagnostix

[www.3ddx.com](http://www.3ddx.com)

### D3 Solutions

[www.d3implantsolutions.com](http://www.d3implantsolutions.com)

### Pro-Art Dental Laboratory

[www.pro-artdentallab.com](http://www.pro-artdentallab.com)

## I Use CEREC

- Set up your software to produce the guide with a hole that matches the Versah® C-Guide™ Sleeve outer-diameter.
- Create buccal window to match the Versah® C-Guide™ Sleeve opening.
- Order the C-Guide™ Sleeve and the compatible G-Stop™ from Versah®.

## I Print My Own Guide

- Set up your software to produce the guide with a hole that matches the Versah® C-Guide™ Sleeve outer-diameter.
- Create buccal window to match the Versah® C-Guide™ Sleeve opening.
- Order the C-Guide™ Sleeve and the compatible G-Stop™ from Versah®.

## I'd Like to Use My Lab

- Versah® is committed to certify your dental lab to produce the Universal C-Guide™ Surgical Guide. Please call us at: 517-796-3932.



## Versah® Guided Surgery System Compatible Software

A full list of compatible software available at [versah.com/guided-surgery-system/](http://versah.com/guided-surgery-system/)



[3shape.com](http://3shape.com)



[blueskybio.com](http://blueskybio.com)



[codiagnostix.com](http://codiagnostix.com)



[planmeca.com](http://planmeca.com)



powered by swissmeda

[mysmop.com](http://mysmop.com)

# Sizing the Versah<sup>®</sup> Guided Surgery System

G-Stop<sup>™</sup> Gauge




Available  
Sizes

**S**

Placement up  
to 3.25mm  
Implant

Densah<sup>®</sup> Bur  
Compatibility



G-Stop<sup>™</sup> Vertical  
Gauge + Key



OD 4.2mm

C-Guide<sup>™</sup>



**M**

Placement up  
to 4.3mm  
Implant



OD 5mm



**L**

Placement up  
to 5.3mm  
Implant



OD 6mm



**XL**

Placement up  
to 6.2mm  
Implant



OD 7mm



C-Guide<sup>™</sup>  
Sleeve



## Available G-Stop<sup>™</sup> Drilling Depths



[www.versah.com](http://www.versah.com)

**Versah<sup>®</sup>**  
THE OSSEODENSIFICATION COMPANY

©2018 Huwais IP Holdings LLC. All rights reserved. Versah<sup>®</sup> and Densah<sup>®</sup> are registered trademarks of Huwais IP Holdings LLC. HT017 REV07