

Vista-Tec LupOn

Protective shields for Loupes

Vista-Tec LupOn kit protects lenses and loupes systems thanks to its wide and durable polycarbonate shield. The wide polycarbonate shields, available in 2 sizes (medium and large), fit into the support clips attached to the frame by the silicone bands oRings. Vista-Tec LupOn clips can also be used on spectacles with temple arms over 1 cm high.

Vista-Tec LupOn selection guide

use Lupon
M-Shield with
Galilean Loupes



use Lupon
L-Shield with
prismatic Loupes



- protect your lenses and loupes system
- premium optical quality shield
- disinfectable, durable and replaceable shield
- autoclavable LupOn clips and oRings

Vista-Tec LupOn kit

REF 5633
1 pair LupOn clips, 1 M-Shield
1 L-Shield, 3 pairs oRings S,M,L

Refill

Vista-Tec lupOn M-Shield
REF 5635 - 5 M-Shields

Vista-Tec lupOn L-Shield
REF 5636 - 5 L-Shields

Vista-Tec LupOn oRings
REF 5634 - 3 pairs oRings S,M,L



Vista-Tec ilGrande

The extended protection

Vista-Tec ilGrande is much more than an elegant face protection equipment. It is a statement of respect towards yourself, your practice and your patients. Its wide, resistant and replaceable polycarbonate shield grants a high level of protection against liquid splashes, while the premium optical properties provide pleasant and relaxed vision. Vista-Tec ilGrande also features a head band which makes it comfortable to wear on longer working days. Both frame and head band can be autoclaved for maximum safety.

Vista-Tec ilGrande kit

REF 5663
1 frame, 2 ilGrande shields, 1 head band

Vista-Tec ilGrande shield

REF 5666
5 polycarbonate shields



- XXL shields for increased protection
- premium optical quality shields
- CE EN 166:2001 certified
- ultra-light design
- disinfectable durable and replaceable shields
- autoclavable frame and head band
- swiss-made





Vista-Tec laStoria

The original face shield
Since 1989

The elegant Vista-Tec face shield, appreciated by dental and health professionals since its conception in 1989, now redesigned with enhanced protection features. The ultra-lightweight face shield includes a durable, replaceable polycarbonate shield that provides premium optical properties for a relaxed vision, reducing unnecessary eyestrain over longer work hours. Vista-Tec laStoria features a silicone head band that enables optimized stability and comfort.

Vista-Tec laStoria kit
REF 5621
1 frame, 2 laStoria shields, 1 head band


Vista-Tec laStoria shield
REF 5626
5 polycarbonate shields

- premium optical quality shields
- disinfectable, durable and replaceable shields
- comfortable head band for ideal stability
- certified CE EN166:2001
- ultra-light, elegant & timeless design
- made in Switzerland



more info at
polydentia.ch

 **polydentia**
swiss manufacture



For spectacles or loupes
with temple arms of max
height of 1cm

Vista-Tec ClipOn

Protective shields for thin-frame
spectacles and loupes




ClipOn enables fixing a protective shield on spectacles. The practical ClipOn clips provide support for both ClipOn shield and the larger L-ClipOn shield, ideal for Galilean loupes. ClipOn clips can be inserted on spectacles or loupes with temple arms of a maximum height of 1cm.

- protect your lenses and loupes system
- ideal for thin-frame spectacles
- disinfectable, durable and replaceable shields

Vista-Tec ClipOn kit
REF 5622
1 pair ClipOn clips, 5 ClipOn shields

Vista-Tec ClipOn shield
REF 5646
10 shields



Vista-Tec L-ClipOn kit
REF 5630
1 pair ClipOn clips, 5 L-ClipOn shields

Vista-Tec L-ClipOn shield
REF 5632
10 shields